


## EnzyChrom™ Ketone Body Assay Kit (EKBD-100)

### Quantitative Colorimetric Determination of Ketone Body at 340nm

#### DESCRIPTION

**KETONE BODIES** (acetoacetic acid and 3-hydroxybutyric acid) are produced in the liver mainly from oxidation of fatty acids, and are normally present at low concentrations in urine and blood. Increased ketone concentrations in the blood may lead to metabolic acidosis, which has been associated with diabetes, childhood hypo-glycaemia, growth hormone deficiency, alcohol or salicylate intoxication and inborn errors of metabolism.

Simple, direct and automation-ready procedures for measuring acetoacetic acid (AcAc) and 3-hydroxybutyric acid (BOH) are very desirable. BioAssay Systems' EnzyChrom™ ketone body assay is based on 3-hydroxybutyrate dehydrogenase catalyzed reactions, in which the change in NADH absorbance, measured at 340nm, is directly related to the AcAc and BOH concentrations,



#### APPLICATIONS

Direct assays of ketone body in serum, plasma, urine and other biological samples.

#### KEY FEATURES

**Sensitive and accurate.** Uses 10  $\mu\text{L}$  sample. Linear detection range of 0.12 to 8 mM for each ketone body in 96-well plate assay.

**Convenient.** The procedure involves adding a single working reagent, and reading the optical density at room temperature.

**High-throughput.** Can be automated as a high-throughput 96-well plate assay for many samples per day.

#### KIT CONTENTS (200 TESTS IN 96-WELL PLATES)

<b>AcAc Buffer:</b>	20 mL	<b>BOH Buffer:</b>	20 mL
<b>AcAc Reagent:</b>	Dried	<b>BOH Reagent:</b>	1 mL
<b>AcAc Standard:</b>	200 $\mu\text{L}$	<b>BOH Standard:</b>	200 $\mu\text{L}$
<b>HBDH Enzyme:</b>	120 $\mu\text{L}$		

**Storage conditions.** The kit is shipped on ice. Store all reagents at  $-20^\circ\text{C}$ . Shelf life: 6 months after receipt, 3 weeks after reconstitution.

**Precautions:** reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

#### PROCEDURES

**Samples:** serum and plasma samples should be non-hemolyzed and assayed immediately. If not assayed, samples can be stored at  $-80^\circ\text{C}$  for up to 30 days.

**Reagent preparation:** bring all reagents to room temperature prior to assay. Reconstitute the AcAc Reagent tube with 1000  $\mu\text{L}$  dH<sub>2</sub>O (final 10 mM). Unused AcAc Reagent is stable for three weeks when stored frozen at  $-20^\circ\text{C}$ . During experiment, keep the HBDH enzyme on ice or in refrigerator ( $2-8^\circ\text{C}$ ).

#### AcAc Assay

1. **Standards.** Prepare 8 mM standard by mixing 5  $\mu\text{L}$  AcAc standard with 45  $\mu\text{L}$  distilled H<sub>2</sub>O. Transfer 5  $\mu\text{L}$  H<sub>2</sub>O, 5  $\mu\text{L}$  8 mM AcAc standard in separate wells of a clear, flat-bottom, 96-well plate.

**Samples.** Transfer 5  $\mu\text{L}$  sample into two wells, one *Sample* well and one sample *Blank* well.

2. **Reaction.** Prepare Working Reagent for H<sub>2</sub>O, Standard and *Sample* wells, by mixing 195  $\mu\text{L}$  AcAc Buffer, 8  $\mu\text{L}$  reconstituted AcAc Reagent and 0.5  $\mu\text{L}$  HBDH Enzyme for each well. The Blank Reagent is prepared by mixing, for each *blank* well, 195  $\mu\text{L}$  AcAc Buffer and 8  $\mu\text{L}$  reconstituted AcAc Reagent (i.e., *no enzyme*).

Add 195  $\mu\text{L}$  Working Reagent to the H<sub>2</sub>O, Standard and *Sample* wells. Add 195  $\mu\text{L}$  Blank Reagent to *Sample Blank* wells. Gently tap plate to mix.

3. Incubate 5 min at room temperature. Read OD<sub>340nm</sub>. Calculate the acetoacetic acid (AcAc) concentration from the OD values of the H<sub>2</sub>O, 8 mM Standard, Sample and Sample Blank wells,

$$[\text{AcAc}] = \frac{\text{OD}_{\text{BLANK}} - \text{OD}_{\text{SAMPLE}}}{\text{OD}_{\text{H}_2\text{O}} - \text{OD}_{\text{STANDARD}}} \times 8 \text{ (mM)}$$

#### BOH Assay

1. **Standards.** Prepare 8 mM standard by mixing 5  $\mu\text{L}$  BOH standard with 45  $\mu\text{L}$  distilled H<sub>2</sub>O. Transfer 5  $\mu\text{L}$  H<sub>2</sub>O, 5  $\mu\text{L}$  8 mM BOH standard in separate wells of a clear, flat-bottom, 96-well plate.

**Samples.** Transfer 5  $\mu\text{L}$  sample into two wells, one *Sample* well and one sample *Blank* well.

2. **Reaction.** Prepare Working Reagent for H<sub>2</sub>O, Standard and *Sample* wells, by mixing 195  $\mu\text{L}$  BOH Buffer, 8  $\mu\text{L}$  BOH Reagent and 0.5  $\mu\text{L}$  HBDH Enzyme for each well. The Blank Reagent is prepared by mixing, for each *blank* well, 195  $\mu\text{L}$  BOH Buffer and 8  $\mu\text{L}$  BOH Reagent (i.e., *no enzyme*).

Add 195  $\mu\text{L}$  Working Reagent to the H<sub>2</sub>O, Standard and *Sample* wells. Add 195  $\mu\text{L}$  Blank Reagent to *Sample Blank* wells. Gently tap plate to mix.

3. Incubate 15 min at room temperature and read OD<sub>340nm</sub>. Calculate the 3-hydroxybutyric acid (BOH) concentration from the OD values of the sample, sample blank, Standard and H<sub>2</sub>O,

$$[\text{BOH}] = \frac{\text{OD}_{\text{SAMPLE}} - \text{OD}_{\text{BLANK}}}{\text{OD}_{\text{STANDARD}} - \text{OD}_{\text{H}_2\text{O}}} \times 8 \text{ (mM)}$$


Total ketone body (TKB) concentration is calculated as,

$$[\text{TKB}] = [\text{AcAc}] + [\text{BOH}]$$

**Note:** if the calculated [AcAc] or [BOH] is higher than 8 mM, dilute sample in H<sub>2</sub>O and repeat this assay. Multiply the results by the dilution factor.

#### MATERIALS REQUIRED, BUT NOT PROVIDED

Pipeting (multi-channel) devices. Clear flat-bottom 96-well plates (e.g. Corning Costar) and plate reader.



Standard Curves of Acetoacetic Acid (AcAc) and 3-Hydroxybutyric Acid (BOH)

#### LITERATURE

- Nuwayhid, N.F., Johnson, G.F. and Feld, R.D. (1988). Kinetic measurement of the combined concentrations of acetoacetate and  $\beta$ -3-hydroxybutyrate in serum. *Clin. Chem.* 34/9, 1790-1793.
- Hansen, J.L. and Freier, E.F. (1978). Direct assays of lactate, pyruvate,  $\beta$ -3-hydroxybutyrate, and acetoacetate with a centrifugal analyzer. *Clin. Chem.* 24/3, 475-479.
- Siegel, L., Robin, N.I. and McDonald, L.J. (1977). New approach to determination of total ketone bodies in serum. *Clin. Chem.* 23/1, 46-49.